Modèle 15007

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| EMC Declaration                                                                                                                                                                       |                                                                                                                                                                                           |                            |                                                                                                                                                                                                                                                                                   |
|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Guidance and manufacturer's declaration – electromagnetic immunity                                                                                                                    |                                                                                                                                                                                           |                            |                                                                                                                                                                                                                                                                                   |
| The "JXB-178" is intended for use in the electromagnetic environment specified below. The customer or the user of the "JXB-178" should ensure that it is used in such an environment. |                                                                                                                                                                                           |                            |                                                                                                                                                                                                                                                                                   |
| Immunity test                                                                                                                                                                         | IEC 60601 test level                                                                                                                                                                      | Compliance level           | Electromagnetic environment – guidance                                                                                                                                                                                                                                            |
| Electrostatic discharge (ESD) IEC 61000-4-2                                                                                                                                           | ±6 kV contact<br>±8 kV air                                                                                                                                                                | ±6 kV contact<br>±8 kV air | Floors should be wood, concrete or ceramic tile. If floors are covered withsynthetic material, the relative humidityshould be at least 30 %.                                                                                                                                      |
| Electrical fast transient/burst IEC 61000-4-4                                                                                                                                         | ±2 kV for power supply lines<br>±1 kV for input/output lines                                                                                                                              | Not Applicable             | Mains power quality should be that of a typical commercial or hospital environment.                                                                                                                                                                                               |
| Surge IEC 61000-4-5                                                                                                                                                                   | ±1 kV line(s) to line(s)<br>±2 kV line(s) to earth                                                                                                                                        | Not Applicable             | Mains power quality should be that of a typical commercial or hospital environment.                                                                                                                                                                                               |
| Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11                                                                                   | <5 % UT<br>(>95 % dip in UT) for 0,5 cycle<br><br>40 % UT<br>(60 % dip in UT) for 5 cycles<br><br>70 % UT<br>(30 % dip in UT) for 25 cycles<br><br><5 % UT<br>(>95 % dip in UT) for 5 sec | Not Applicable             | Mains power quality should be that of a typical commercial or hospital environment. If the user of the "JXB-178" requires continued operation during powermains interruptions, it is recommended that the "JXB-178" be powered from an uninterruptible power supply or a battery. |
| Power frequency (50/60 Hz) magnetic field IEC 61000-4-8                                                                                                                               | 3 A/m                                                                                                                                                                                     | 3 A/m                      | Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial hospital environment.                                                                                                                                            |
| NOTE UT is the a.c. mains voltage prior to application of the test level.                                                                                                             |                                                                                                                                                                                           |                            |                                                                                                                                                                                                                                                                                   |

## **DeltaTrak Corporate**

📍 P.O. Box 398 Pleasanton CA 94566

EO: sales@eofirm.com | www.eofirm.com | 0755-84200058

EO: sales@eofirm.com | www.eofirm.com | 0755-84200058

## **DTI Mexico International**

📍 Guadalajara, Mexico

EO: sales@eofirm.com | www.eofirm.com | 0755-84200058

EO: sales@eofirm.com | www.eofirm.com | 0755-84200058

## **DTI South America SA**

📍 Santiago, Chile

EO: sales@eofirm.com | www.eofirm.com | 0755-84200058

EO: sales@eofirm.com | www.eofirm.com | 0755-84200058

## **DTI China Limited**

📍 Shenzhen China

EO: sales@eofirm.com | www.eofirm.com | 0755-84200058

2031-2004 | 0923-2110

EO: sales@eofirm.com | www.eofirm.com | 0755-84200058

## **DTI Europe bvba**

📍 Antwerp Belgium

EO: sales@eofirm.com | www.eofirm.com | 0755-84200058

EO: sales@eofirm.com | www.eofirm.com | 0755-84200058

## **DTI Japan Limited**

📍 Osaka, Japan

EO: sales@eofirm.com | www.eofirm.com | 0755-84200058

EO: sales@eofirm.com | www.eofirm.com | 0755-84200058

## **DTI Asia Pacific**

📍 NT, Hong Kong

EO: sales@eofirm.com | www.eofirm.com | 0755-84200058

EO: sales@eofirm.com | www.eofirm.com | 0755-84200058

## **DTI South Pacific**

📍 Auckland, New Zealand

EO: sales@eofirm.com | www.eofirm.com | 0755-84200058

EO: sales@eofirm.com | www.eofirm.com | 0755-84200058